
## A multi-compartmental model of the HepQuant SHUNT test quantifies anatomic shunting and stratifies risk for varices: combined results from the HALT-C and SHUNT-V studies <u>M.P. McRae<sup>1</sup>, S.M. Helmke<sup>2</sup>, and G.T. Everson<sup>2</sup></u> <sup>1</sup>Custom Diagnostic Solutions LLC, Houston, TX, USA; <sup>2</sup>HepQuant LLC, Denver, CO, USA Figure 4. Compartmental and Minimal Model DSI Correlation **1. Introduction** 4. Results The compartmental model fit the data with median adjusted R<sup>2</sup> of 0.99 for the IV 13C-CA curves and 0.95 for the Oral d4-CA curves. The disease severity $R^2$ =0.9912 50 Significant differences were observed in averaged oral and IV curves between controls, no varices, index (DSI) calculated Slope=0.9822 small varices, and large varices, in plasma concentration-time profiles (Figure 2). from the SO 40 intercept=-0.8092 compartmental model • These results suggest the HepQuant test is detecting the presence and the size/magnitude of anatomic IV and oral curves shunts and collateral circulation in patients with varices. Θ correlates well to the **Ö** 30 Figure 2. Population-averaged IV curve (left) oral curve (right) simulations **DSI calculated from the** (mean and 95% confidence) from the compartment model. minimal model IV and <u>a</u>20 **2.** Aim oral curves, which was IV Oral 0 0 10 previously validated as 10 Large Varices a measure of global liver function [4] -Small Varices (Figure 4). -No Varices -Control 0 3. Methods Min. Model DSI $(\mathbf{r})$ Con **5.** Conclusions 150 100 150 50 Time (min) Time (min) with validated indices of hepatic disease. Plasma concentrations normalized to 75 kg bodyweight. Inset on IV curve is magnified view of 20-60 minutes. $\succ$ The median anatomic shunt flow rate (q<sub>PS</sub>) demonstrated significant differences between all groups (Wilcoxon rank-sum p < 0.05) (Figure 3, left). HepQuant's global liver function tests. In parallel, the reduced portal inflow to the liver (q<sub>PT</sub>) was observed among varices groups, with significant difference between all groups (p < 0.001) (Figure 3, right) **6.** References Figure 3. Scatter boxplots of compartmental parameters representing anatomic shunt and portal inflow to the liver across patient groups (control, [1] Lee, W.M., et al., Evolution of the HALT-C Trial: pegylated interferon as maintenance none, small, and large esophageal varices). therapy for chronic hepatitis C in previous interferon nonresponders. Controlled Clinical Trials, 2004. 25(5): p. 472-492. **Portal Inflow to Liver Anatomic Shunt** [2] Everson, G.T., et al., The spectrum of hepatic functional impairment in compensated 1.5 chronic hepatitis C: results from the Hepatitis C Anti-viral Long-term Treatment against Oral d4-CA Cirrhosis Trial1. Alimentary Pharmacology & Therapeutics, 2008. 27(9): p. 798-809. [3] The SHUNT-V Study for Varices, ClinicalTrials.gov study ID: NCT03583996. $\mathbf{X}$ \_ [4] Everson, G.T., et al., Portal-systemic shunting in patients with fibrosis or cirrhosis due in flow rates 0.6 to chronic hepatitis C: the minimal model for measuring cholate clearances and shunt. between Alimentary Pharmacology & Therapeutics, 2007. 26(3): p. 401-410 $q_{SP}$ し compartments 0.4 Portal ഗ് വ.5 7. Disclosures $q_{PS}$ MPM is a paid consultant for HepQuant LLC. SMH and GTE are employees and equity $q_{PT}$ Cl<sub>H</sub> members of HepQuant LLC. All authors have provisional patents pending. HepQuant tests are not FDA approved and are for investigational use only under FDA guidelines for total hepatic Liver investigational device exemption (IDE).





In chronic liver disease (CLD), persistent inflammation and fibrosis impair hepatocyte function and increase portal pressure, resulting in portal venous collaterals and portalsystemic shunting. The HepQuant SHUNT test uses stable nonradioactive 24-<sup>13</sup>C-cholate (13C-CA) intravenously and 2,2,4,4-d4-cholate (d4-CA) orally to simultaneously measure clearances from both portal and systemic circulations.

To determine whether compartmental model of the HepQuant SHUNT test could quantify anatomic shunting in CLD patients with esophageal varices.


- > A compartmental model was developed using HepQuant SHUNT test results from the HALT-C Trial [1,2], the SHUNT-V Study [3], and lean controls: in total, 492 CLD patients with esophago-gastro-duodenoscopy findings (291 no varices, 143 small varices, 58 large varices) and 26 lean controls.
- A compartmental model (Figure 1) describes dual IV and oral cholate clearance and is divided into 3 compartments: Systemic, Portal, and Liver.
- Transfer between compartments was modeled by a system of 18 differential equations with assumptions from measured and literature-derived values.
- Model parameters were solved for each subject estimated by nonlinear least-squares regression including the anatomic shunt flow (q<sub>PS</sub>) and portal venous inflow to the liver (q<sub>PT</sub>).
- Simulations and parameter estimation was solved for each patient on the IV and oral cholate plasma concentrationtime data simultaneously.







Portal inflow to liver was normalized to 75 kg bodyweight.



The compartmental model discriminated patients with small and large varices, attributed cholate clearance to hepatocyte function and anatomic shunting, and correlated

Compartmental analysis has significant potential to enhance diagnostic performance and clinical utility of

# 8. Contact Information

Gregory T. Everson: greg.everson@hepquant.com Steve M. Helmke: steve.helmke@hepquant.com Michael P. McRae: mpmcrae@customdxsolutions.com